






## MOTIVATION CHILDREN'S ACTIVE WHEELCHAIR

**MANUFACTURER:** Motivation

**OVERVIEW:** An adjustable, active style wheelchair for children. Designed for better performance over rough ground and in conditions normally found in rural areas, but can also be used in urban and semi-urban area. A single medium size rubber castor wheel at the front and a long wheelbase helps the user to propel safely over uneven and inclined ground. Includes a contoured pressure relief cushion, which can be cut to size for each user.

**INTENDED USERS:** Children without complex postural support needs who regularly use a wheelchair



|   |  |  |
|---|--|--|
| PRODUCT CODE  | K3 - 02  |  |
|  SEAT WIDTH        | 12" (30 cm)  |  |
|  SEAT DEPTH        | 9.4" (24 cm), 10.4" (29 cm), 13.4" (34 cm), adjustable   |  |
|  FOOTREST HEIGHT | Adjustable, 8.7" – 11.4" (22 cm – 29 cm) measured from the seat to the footplate   |  |
| BACKREST HEIGHT   | 7.9" (20 cm), 14" (26.5 cm), 13" (33 cm), adjustable   |  |
| BACK ANGLE  | Fixed  |  |
| ARM SUPPORT   | Fixed  |  |
| REAR WHEEL DIMENSIONS   | 20" x 1.75" diameter (50.8 cm x 4.4 cm) pneumatic tires, with 36 spoke aluminum rims   |  |
| CASTER DIMENSIONS   | 6" x 1.6" diameter (15 cm x 4 cm) rubber castor wheel  |  |
| PRODUCT WEIGHT  | 24.9 lb (11.3 kg) including cushion  |  |
| WEIGHT CAPACITY   | 110 lb (50 kg)   |  |
| OTHER FEATURES  | <ul style="list-style-type: none"> <li>• Rigid frame</li> <li>• Padded fabric backrest with tension-adjustable straps</li> <li>• Adjustable backrest height</li> <li>• Adjustable footrest height, depth, and angle</li> <li>• Adjustable seat depth</li> <li>• Height adjustable push handles</li> <li>• Calf strap</li> <li>• Pelvic strap</li> <li>• Anti-tip wheels</li> <li>• Pneumatic, medium tread tires</li> <li>• Chromium plated steel pushrims</li> <li>• Quick release axles</li> <li>• Includes a contoured pressure relief cushion</li> <li>• Easily collapsible for transport</li> <li>• Certifications: ISO 14001, ISO 13845, OHS 10001, passed ISO 7176, including fire retardant foam and fabric</li> </ul> |  |
| WARRANTY  | All registered wheelchairs receive 12 months on the frame and 6 months on relevant spare parts. See <a href="https://clasp.org/warranty">clasp.org/warranty</a> for additional information.  |  |
| ACCESSORIES   | <ul style="list-style-type: none"> <li>• Tire repair kit</li> <li>• Tire pump</li> <li>• Basic tool kit: two multi-function spanners, one 4 mm Allen key, two 5 mm Allen keys, one 6 mm Allen key</li> </ul>   |  |



9509 Vassar Avenue, Los Angeles, CA 91311, USA  
TEL: +1 818.255.0100

[WWW.CLASPHUB.ORG](http://WWW.CLASPHUB.ORG)

*This product information was produced by Momentum Wheels for Humanity through Consolidating Logistics for Assistive Technology Supply and Provision II (CLASP II), a four-year cooperative agreement funded by the U.S. Agency for International Development under Agreement Number AID-OAA-A-17-00021, beginning August 30th, 2017.*

*CLASP II is a USAID-funded project implemented by Momentum Wheels for Humanity that increases access to a variety of high-quality, affordable assistive products, and promotes appropriate provision globally.*