



## MOTIVATION CLIP-ON TRICYCLE

**MANUFACTURER:** Motivation Charitable Trust


**OVERVIEW:** The Clip-on Tricycle attachment is a quick release clip-on, clip-off attachment that can be fitted to the Motivation Rough Terrain wheelchair. It converts the wheelchair into a tricycle allowing the user to travel greater distances more efficiently.

**INTENDED USERS:** Active adults and teens whose lifestyle includes traveling long distance or going over rough terrain.



PRODUCT CODE	WM3TRCO-01
WHEELS	20" x 1.75" diameter (51cm x 4cm), pneumatic front wheel with 28-hole steel rim
PEDAL POSITION	The pedal position can be adjusted horizontally by approximately 2" (5 cm); vertical adjustment is fixed and is designed to suit most users
PRODUCT WEIGHT	35.3 lbs (16 kg)
BRAKES	Strong V-brake on the front wheel
CHAIN	Standard bicycle chain is included; chain tension is adjustable
GEARS	Simple single speed freewheel gear to suit most users
CERTIFICATIONS	ISO 9000
ADDITIONAL INFORMATION	Includes tricycle attachment <b>only</b> . Does not include Motivation Rough Terrain wheelchair.

**BEST FOR USER ENVIRONMENT:**





9509 Vassar Avenue, Los Angeles, CA 91311, USA  
TEL: +1 818.255.0100

[WWW.CLASPHUB.ORG](http://WWW.CLASPHUB.ORG)

*This product information was produced by Momentum Wheels for Humanity through Consolidating Logistics for Assistive Technology Supply and Provision II (CLASP II), a four-year cooperative agreement funded by the U.S. Agency for International Development under Agreement Number AID-OAA-A-17-00021, beginning August 30th, 2017.*

*CLASP II is a USAID-funded project implemented by Momentum Wheels for Humanity that increases access to a variety of high-quality, affordable assistive products, and promotes appropriate provision globally.*

